

## Carbonyl Unknowns

### Overview:

You will receive a carbonyl compound as an unknown. It can be either an aldehyde or a ketone, and may or may not contain an aromatic ring. Your job will be to identify your carbonyl compound. Several pieces of information will be useful:

- NMR information on the starting material. (H, C-decoupled, and 2D HC-NMR)
- Water solubility tests (big or small? Aromatic or not?)
- Boiling point of starting material (try at least once)
- The melting point of the derivative (required)


### Classifying Tests

1. Water Solubility Test (Helpful, but not always decisive or clear-cut. Use, but don't depend on it too much?!)
  - Add 15 drops of water to a small test tube, and then add 2 drops of sample. Stir vigorously. Is it homogeneous or heterogeneous? If heterogeneous, do the droplets float or sink?
  - Interpretation: **Insoluble proves  $\geq 3$  carbons. Soluble proves  $\leq 6$  carbons.**
    - a. **Small:** Carbonyls with  $< 4$  carbons always **dissolve**
    - b. **Big:** Carbonyls with  $> 6$  carbons **never dissolve**
    - c. **Borderline:** Carbonyls with **4-6 C's, borderline, may dissolve or may not.** Sometimes adding some more water will dissolve, if doesn't initially.
    - d. An **insoluble carbonyl that sinks** has an **aromatic** ring present for sure
    - e. An **insoluble carbonyl that floats** is **probably nonaromatic**, although some aromatics are also floaters.
2. Summary of chemical tests related to carbonyls, (we'll only do 2,4-DNP), but which you should know to answer post-lab questions
  - 2,4-dinitrophenylhydrazine (**DNP**) test: positive for **aldehydes or ketones**.
  - **Tollens'** test: Positive for **Aldehydes**, not for Ketones. Similar to Schiff's test, but more famous (good) but more expensive (bad) (Note: we'll tell by H-NMR. Shift at  $\sim 9$ -10ppm proves aldehyde.)
  - **Iodoform** Test: Positive for **Methyl Ketones** ( $\text{CH}_3\text{COR}$ ). (This is also pretty easy to see by H-NMR, since you get a 3H singlet in the 2's.)
  - **$\text{Br}_2/\text{CH}_2\text{Cl}_2$**  test: Positive for **Alkenes** (to distinguish  $\text{C}=\text{C}$  from  $\text{C}=\text{O}$  double bonds)
3. **NMR: Prepare a sample** by drawing up about 1 inch of your unknown into the skinny part of a long-stemmed pipet, then place the pipet into an NMR tube. Add 0.8 mL of  $\text{CDCl}_3$  directly through the pipet to rinse the sample into the NMR tube. Cap it, shake it, and set it into the NMR queue and run the experiment called "H\_C\_HC" in the 355-365 folder.
  - **Aldehydes** are easily distinguished from ketones by H-NMR. The aldehyde hydrogen, which is attached to the carbonyl carbon, **shows up in the 9-10.5ppm chemical shift area.** Ketones will show no such signal in that area.
  - **Aromatic hydrogens ortho to a carbonyl are typically pushed downfield, toward 8 ppm.** This is because a carbonyl group is a strong electron withdrawer, so it makes the ortho carbons more electron poor, which "deshields" the ortho hydrogens.




**SOME CHEMICAL TESTS TO KNOW**


2,4-Dinitrophenylhydrazine (“DNP”) Test: Specific for Aldehydes or Ketones (but not esters, acids, or amides)



The “DNP” test is positive for both aldehydes and ketones, but not for alkenes or esters/acids/amides. This is representative of how  $\text{H}_2\text{N}-\text{Z}$  reagents react with aldehydes or ketones to eliminate water and make “imines”, with a  $\text{C}=\text{N}-\text{Z}$  bond. In the chemical test, the DNP reagent is soluble; if a derivative forms, it precipitates from solution. So the formation of a precipitate is what you watch for. The DNP-derivatives tend to be highly crystalline because of the extended conjugation; from the carbonyl carbon through the two nitrogens through the ring through the two nitro attachments, all the atoms are flat and  $\text{sp}^2$ . The color of the precipitate is often informative; saturated carbonyl compounds tend to give yellow derivatives, while unsaturated aldehydes or ketones tend to give red or orange derivatives. The experiment is excellent as a chemical test, when you don’t know if you have an aldehyde or ketone. But it is also excellent as a way to make a solid derivative which can be purified by recrystallization and whose melting point can be taken. The melting points of many DNP derivatives are known and listed.

**Tollens Test: Specific for Aldehydes. Positive for Aldehydes Only.**


A classic alternative to the Schiff’s test for aldehydes is the Tollens Test. Tollens reagent is a soluble  $\text{AgOH}$  solution. [Actually  $\text{Ag}(\text{NH}_3)_2\text{OH}$ ]. When mixed with an aldehyde, the aldehyde carbon is oxidized to a carboxylic acid, and the  $\text{Ag}(\text{I})$  cation is reduced to elemental  $\text{Ag}(0)$ . The elemental silver films out on the surface of the test tube in which the test is conducted, and a “silver mirror” can be observed. This reaction has historic importance. For centuries during the middle ages this was the process used to make mirrors. (These silver mirrors were less clear than modern mirrors). This silver coating process was also used to apply a silver coating to any object. We will not use this test in lab because the Schiff’s test is cheaper and easier. Test tubes used for Tollens’ test must be thrown away, and the silver reagent is somewhat expensive.

**Iodoform Test: Specific for Methyl Ketones (CH<sub>3</sub>COR)**

Methyl ketones can be distinguished from other ketones by the iodoform test. The methyl ketone is treated with iodine in an NaOH/water solution. Methyl ketones produce a yellow solid called “iodoform”, other ketones or aldehydes do not. The mechanism is shown below, and is somewhat complex. Deprotonation of the methyl ketone hydrogen gives a resonance-stabilized anion, which attacks iodine. Once the first iodine is installed, the remaining methyl hydrogens become even more acidic and get deprotonated followed by iodination in rapid sequence to generate the tri-iodo RCOI<sub>3</sub> species (in box). Hydroxide routinely adds to carbonyls, but normally this addition is reversible, non-productive, and insignificant. However, hydroxide addition to the RCOI<sub>3</sub> is productive; in this case, the anion (in circle) can eliminate the <sup>-</sup>Cl<sub>3</sub> anion. This is a decent leaving group because the three electron-withdrawing iodo groups stabilize the anion. This elimination is also irreversible, so by LeChatelier’s principle all of the chemicals drain off through this pathway. Following elimination, the <sup>-</sup>Cl<sub>3</sub> anion picks up a proton to make iodoform, CHI<sub>3</sub>, which is a yellow crystalline solid. The formation of this yellow solid is a “positive” test; if no yellow solid forms, the test is “negative”. Ketones other than methyl ketones are unable to get to the RCOI<sub>3</sub> species (in box), are unable to undergo the fragmentation that the circled anion undergoes, and are unable to make the solid iodoform.

**Br<sub>2</sub> Test: Specific for Alkenes (Not Ketones or Aldehydes)**

Bromine is a routine test for alkenes. (Although a mono-substituted alkene is shown in the picture, di-, tri- and tetra-substituted alkenes also react with bromine.) Bromine adds to alkenes but not to carbonyl compounds (or to ordinary arenes). The nature of the test



is to add a few drops of bromine, which is strongly colored, to an excess of an organic sample. If the color disappears, it means the bromine reacted and therefore that the organic unknown contains an alkene. If the color persists, it means the bromine did not react, and therefore that no alkene is present in the organic unknown.

Aldehyde/Ketone Candidates

Bp of Starting Carbonyl	Unknown	mp of 2,4-DNP Derivative
48	propanal	148
56	acetone	126
63	2-methylpropanal	187(183)
75	butanal	123
80	2-butanone	117
91	3-methylbutanal	123
92	2-methylbutanal	120
100	2-pentanone	143
102	3-pentanone	156
103	pentanal	107(98)
115	4-methyl-2-pentanone	95
128	5-hexen-2-one	108
129	4-methyl-3-penten-2-one	205
131	cyclopentanone	146
131	hexanal	104(107)
145	4-heptanone	75
145	5-methyl-2-hexanone	95
146	2-heptanone	63-68 *
147	3-heptanone	81
153	heptanal	108
156	cyclohexanone	162
169	3-methylcyclohexanone	155
173	2-octanone	58
179	benzaldehyde (PhCHO)	237
200	o-methylbenzaldehyde	194
204	p-methylbenzaldehyde	234
202	ethanoylbenzene	244
216	1-phenyl-2-propanone	156
217	(2-methylpropanoyl)benzene	163
218	propanoylbenzene	191-198 *
226	p-methylacetophenone	258
232	butanoylbenzene	191
235	4-phenyl-2-butanone	127
248	p-methoxybenzaldehyde	253



**Name:**

**Lab Report Requirements:** No procedure or yield information required. Fill out the unknown report sheet. Attach your NMRs. (Must take at least one of H-NMR or C-NMR, or both.) Answer the following questions.

**Questions:**

1. What is the purpose of making derivatives of liquid unknowns?
2. Using a chemical test or tests, how could you distinguish between 3-pentanone and pentanal?
3. Using a chemical test or tests, how could you distinguish between 3-pentanone and 2-pentanone?
4. Using a chemical test or tests, how could you distinguish between 3-pentanone and 4-penten-1-ol?
5. Draw a possible structure for a molecule  $C_5H_8O$  that gives a positive tollens' test and does not react with  $Br_2/CH_2Cl_2$ ?
6. Draw the structure of a compound  $C_5H_8O$  that reacts with 2,4-dinitrophenylhydrazine, decolorizes bromine in dichloromethane, but does not give a positive iodoform test.
7. Draw two structural isomers for  $C_5H_{10}O$  that would both give positive iodoform tests?
8. Draw a possible structure for  $C_4H_8O$  that would not give a positive dinitrophenylhydrazone test?

## Unknown Report Sheet-Carbonyls

Name

Your unknown Letter/Number:

Draw the structure for your unknown:

## 1. Solubility Tests on Starting Material

Solubility in Water: \_\_\_\_\_ If Insoluble, Does it Float or Sink?

Conclusion:

## 2. Boiling point:

3. Derivative:                      observed mp                      literature mp (see p41)

Crude (if possible):

Recrystallized

## 4. H-NMR (attach, with assignments/interpretation.)

- **On the proton spectrum, create a STANDARD 4-column SUMMARY REPORT of your ACTUAL H-NMR data, detailing chemical shifts, integrations, and splittings, and "source".**
- Chemical shifts need to be specified to at least the nearest 0.1 ppm
- Draw the structure of your molecule, with identifiers by each carbon (a, b, c...).
- Then on your standard summary table add a "source" column in which you explain which hydrogens (a, b, or c, etc...) are responsible for which signals. Note: if the sample is too concentrated, the splitting may in some cases get broadened and become problematic.

## 5. C-NMR (attach, with assignments/interpretation)

- Draw the structure of your molecule, with identifiers by each carbon (a, b, c...).
- Draw letters next to carbon lines that can be assigned.
  - The C-H 2D NMR will help you to be able to assign all of your lines.)

## 6. Comments, difficulties, complaints, etc.