Scheme 2 Part F: NMR and GC-MS.

- Checking your NMR prior to starting Scheme 3 is wise.
- Each NMR and each GC/MS will again need to have a clear label that corresponds to the identity of the chemical.
- You MUST again draw actual structures of your product on both NMRs and GC/MS's.
- When discussing NMR's or GC/MS's in the reports, always do so by label. (For example, "GC-MS-6c was taken", or "as shown in NMR-6C, the product was very pure....")

<u>NMR Analysis/Interpretation</u>. Concepts and expectations are similar to the Scheme 1 report.

- 41. An <u>abbreviated summary report</u> process will again be required. Draw the structure and label the different carbons. Then make a table with the chemical shifts for the actual <u>non-aromatic</u> <u>C-H's</u>, and by each one write the letter of the carbon to which it is attached. This will demonstrate that you have analyzed and understand your spectrum.
 - Also include integration.
 - But you do not need to analyze/report the splitting (although you may do so.)
 - Which signal is from the β -H, which are from the α -H's, and which are the benzyl H's?
 - Which signals are the methyls or methoxy signals in **3b**, **3c**, and **3e**?
 - Does your product look reasonably pure?

42. Does your NMR confirm that the reaction, and the solvent-removal basically worked?

- Does it look like your starting reactant **3** is still present, or gone?
- Does it look like there is significant amount of solvent left? (Dichloromethane leaves a singlet at about 5.3 ppm.)
- 43. <u>Chirality, H-non-equivalency, and chemical shift</u>: The chirality of the β-carbon not only makes the two β-H's non-equivalent, but also makes the two hydrogens on the newly attached CH2 carbon nonequivalent. The two benzyl hydrogens should each appear as two doublets.
- 44. <u>Chirality, H-non-equivalency, and splitting</u>: As in product 3, the non-equivalence of the two α -hydrogens, and now the two benzyl hydrogens, complicates their splitting.
 - Each of the two α -hydrogens will usually appear as a four-line "doublet of doublets".
 - The β -H will typically also look like a 4-line "doublet of doublets", unless further split by the methyl group in **3e**.
 - The two benzyl H's are each split by each other, so each should look like a doublet.
 - The appearance of these two new doublets is very diagnostic for product 6 formation!
- 45. For your reports, <u>account for all and only the hydrogens connected to sp³ carbons</u>. You don't need to discuss/present N-H hydrogens or aromatic hydrogens. (<u>There are so many</u> overlapping aromatic H's that they aren't interpretively useful in this case.)
- 46. **Signature signals**: All of the samples will have the interesting α and β -hydrogens, and the benzyl hydrogens (5 hydrogens combined). But there will be other additional signature methyl peaks for **6c**; **6b**; and **6e**.
- 47. Chemical shift logic:
 - The α -hydrogens, being next to a carbonyl, but being also β to nitrogen and perhaps also an aromatic, should fall in the high 2's or perhaps the low 3's.
 - The β -hydrogen in **3a-3d** should show up around the low 4's. For the β -hydrogen in **3e**, we'd anticipate the β -hydrogen to show up in the 3's.
 - For the benzyl hydrogens, they are on a carbon that has both a nitrogen (+2) and a benzene (+1) attached, so we'd expect them to come around the low 4's or high 3's as well.
- 48. **Impurities/contamination**: Recognizing which signals come from the desired product and which do not is again significant.
 - Remember that there should be a logical integration ratio for the main H's in product 6

- 49. **Comparison to Other NMR's**: It may be very interesting to look at how your NMR **6** looks compared to how other NMR's look.
 - How different is your **6** from the **3** that you began with in Scheme 2?
 - How different is your **6** compared to classmates who made different versions of **6**?
 - How clean is your NMR compared to that of classmates who made the same version of 6?
- 50. <u>GC-MS:</u> Clearly label each page of each GC/MS printout.
- 51. Draw the structure and molecular weight for your specific product on each GC-MS sheet.
- 18. **Retention time**? What is the retention time for your **6**? How much longer is it than **3**?
- 19. **Purity:** How pure is your **3** by GC?
- 20. **Mass Spec and Molecular Ion**: For your major product **3**, check in the mass-spec whether there is a molecular ion peak that matches the molecular weight for your product.
- 21. **Mass Spec and Fragmentation**: The weakest break-point is at the N-benzyl bond. You should be able to see a benzyl fragment (PhCH2+ = 99) and a fragment that it molecular weight 99. Do you see both of those fragments?
- 22. **Lab report**: In your lab report, make sure that you have not only attached the labeled GC-MS information, but that you also discuss/present the retention time and purity.

Scheme 2/Week 2 Lab Report:

- 1. Write a standard synthesis style lab report for your Scheme 1 reaction $(3 \rightarrow 6)$.
- 2. Make sure that all structures are drawn explicitly.
 - As always for a synthesis style report, you'll want to draw out the reactants and the products. In this case, be sure you draw your **actual** reactant **3** and product in your reaction.
 - None of your pictures should have an "R1": you should illustrate each structure with your actual R1 group drawn, whether that's methyl or phenyl or 4-methoxyphenyl or whatever.
- 3. Show all calculations. (Including any mole => mass for reactants or products)
- 4. When listing your chemicals/reactants and showing gram \rightarrow mol calculations, make sure that you include your main reactant 3!
- 5. Include procedural details and observations as usual.
- 6. Calculate mass yields, and percent yields, etc., for product 6.
- 7. Include your NMR-6 and GC-MS-6, with clear labels, structures drawn, and the abbreviated summary report for non-aromatic C-H hydrogens for NMR-6.
- 8. Print and attach mass spectra for GC-MS-6.
- 9. Include a results/data/discussion/analysis section. The analysis/discussion section needs to address what the yield information told you, and what the NMR and GC-MS data tells you about both the success and the efficiency of your reaction, and the purity of your product **6**.
- 10. The results/data/discussion/analysis section should summarize what the mass/yield/NMR/GC-MS data is, and what conclusions can be drawn from them. Just attaching the NMR's and GC-MS's without discussing or showing that you understand them will not be good. What is the summary for the key non-aromatic C-H hydrogens in your NMR? What is your GC-retention time? Between the NMR and the GC, did it look like the react 3 was successfully converted to product 6, and does your product 6 look reasonably clean? Or is it obviously significantly contaminated? Was the yield respectable, or terrible?
- 11. Note: Keep extra copies of your NMR and your GC-MS's. Pyrazolidinone **6** functions as the product in week two report, but it is the reactant in the Schemd 3/Week 3 report. So when writing up and analyzing Scheme 3, you'll need information about mass, molecular weight, structure, and mmol of your reactant **6**. You'll also need to have NMR and GC for **6** so that you'll be able to compare your product **10** to reactant **6** and tell whether the reaction really worked. You'll also want copies of **6** for your Final Report after week 3.