

Quick Summary: 1. Sample into Autosampler (1-3) 2. Select Experiment (#6) 3. Submit (#7) (from submit mode)
4. Process + Plot (from spectrometer mode)

User's Guide to NMR: General

- For help, see Dr. Jasperse, Hagen 407J, phone 477-2230

Draft 10/28/15

Default conditions:

- Assumes CDCl₃ as Solvent.
- Want ≥0.8 mL of CDCl₃ normally
- Printing on MSUM printer
- H-NMR

1. Add sample to a Spinner/Turbine

2. Adjust depth by placing the turbine into the golden depth finder

3. Load sample/turbine into autosampler.

- Press the round white **Access Request Button** on the panel below the sample trays/doors
- Wait until "status" light turns to a solid yellow, and the message panel reads "door unlocked"**

4. Opening Program on Computer: Usually already open, and usually to correct "operator"

- If not open: Operator: your class or research group Password: none.
- To switch operator, click **Logout** (below main display screen) and select the correct operator

5. "Submit" vs "Spectrometer" modes: New Study/Submit Queue to submit; Spectrometer to print/view

- Spectrometer => New Study: click "**New Study**" button (lower left)
- Submit queue => Spectrometer: click "**Cancel**" button (lower left)

6. Experiment Selection (from within Submit mode). Usually preselected for organic labs.

- If not already in New Sample/submit queue mode, push **New Study** button on lower left
- Proton8** is the normal H-NMR experiment, under the "UserStudies" folder
- For some classes/operators, **Proton8** has been set to open by default, since most NMR's are regular H-NMR's
- Add experiments as needed from the Experiment Selector.
- To **edit or delete**: **right click** on experiment and select "Open Experiment" or "Delete Experiment"

7. 3 Step Submission (assuming the experiment already specified, and still/already in Submit mode).

- Fill Sample Name** (for both computer filing and printout recognition)
 - Click Sample Spot**: Click on the button showing your sample site. (Remember/record! ☺)
 - Submit**: clicking the **red Submit** button on the lower left side.
- Note**: Can repeat this 3-step sequence for new samples/new students, if running same experiment
 - Comment box**: (can add comments for the paper printout). (Control C to cut and Control K to paste)
 - Other submission options** of possible use for advanced labs, research, or offsite Concordia users:
Solvent; **Autoplot** (offsite Concordia users should turn this off); **Email**; **Email Address** (offsite Concordia users should set this correctly!☺); **Lock**: (with non-deuterated solvent run unlocked), **Shim** (with non-deuterated solvent run 1H PFG); **Tune**

8. Opening Completed Samples for Printing and Processing. ("Spectrometer Mode" required)

- Must be in "Spectrometer" mode, not "Submit" mode. (If already in spectrometer mode, "Study Queue" on left will either display "Location queue" or else "Spectrometer")
 - If in "Submit" mode, "Study Queue" on left will display "Submit queue". Click "**Cancel**" to toggle into "Spectrometer" mode
- In "Spectrometer" mode, must **have "Zones" map displayed** (96 sample nodes show). If not, click on little circle icon (⊞) to the upper left of the spectra-display panel.
 - Right click on sample number**
 - Click "Show Study"
 - Click on file folder name located **on the left**
 - Then **double click on spectrum** you want to view to load it into the spectra-display viewscreen.
 - Process > **Auto Plot or Print**. See next page for more detailed printing and processing instructions.
- Re-click the little circle icon (⊞) to get back to zone map in order to open other files
- To return to "Submission" mode in order to run more samples, click "New study"

9. Logout: Click "Logout" button underneath spectrum-display




Concordia Offsite:

- Print Preview
- Enter Correct Email address


10. Plotting (when wanting non-automatic plots)

- Must be in the process mode. (Highlight “Process” beneath the spectrum display)
- a. Click “Auto Plot” or “Print” button, way on lower right corner of page.
 - b. Re-click if you want to print additional copies for the other students
 - Note to offline Concordia users: this “plot” command will print to MSUM NMR-room printer. ☺
 - For advanced labs or research groups, additional plot preferences are available in the process mode by clicking “Plot (Beneath spectrum display, 2nd from bottom underneath “Start”)

11. Horizontal Expansions

- With spectrum displayed on screen, use a panel of display icons on the far right.
- a. Click on the magnifying glass icon (6th icon down, )
 - b. Move your cursor to the left end of the zone you want to expand, then **hold down left mouse button** and slide it to the other end of the zone you want to expand.
 - To return to the full display, you can either click on the 3rd icon () or the 5th icon ()
 - If the lines aren’t tall enough, type “vsadj” (vertical scale adjust) on the command line.

12. Manual Integration: Defining Integrals Yourself (see #13 to also give nice integral numbers)

- With spectrum displayed, **must be in the process mode** (“Process” beneath the spectrum display)
- a. Choose “Integration” (Beneath spectrum display towards left, 2nd underneath “Start”)
 - b. Hit “Clear Integrals” button (slightly further to the right and lower down from previous button)
 - c. Hit “Interactive Resets” button (immediately above the “clear integrals” button) and define
 1. Move cursor beyond the left end of the signal you want to integrate.
 2. **Left-mouse click-and-release**
 3. Move the cursor to the right of the signal, and again click-and-release. Everything between the two “clicks” will be integrated.
 4. Repeat this for each area you want to integrate.
 - d. Click very top cursor icon () to the right of the display screen to regain normal cursor function

13. Setting Nice Integral Numbers (While already in integration mode following steps a-d above)

- a. Click cursor on one of your integral regions
- b. Click “Normalize Area to” “Single Peak” below “Set Integral Area” panel underneath the display
- c. Set “integral area” to some nice whole number (1, 2, or 3, depending on your molecule)
- d. Click the “set integral value” button
 - If it says “cursor is outside of integral region”, then reset the cursor on an integral of choice, and re-click the “set integral value” button again.
 - Click “Auto Plot” (lower right) in order to print.


14. Other Processing Options for Advanced Users/Research Groups/2D-NMR

- | | | |
|---------------------|---------------------|--|
| 1. Peak Picking | 4. Insets | 6. Absolute Concentration Integration |
| 2. Vsadj | 5. Arraying spectra | 7. 2D NMR processing, including varying the signal intensity |
| 3. wp=2p sp=2p plot | | |

15. Opening Spectra From the Data Folders

- Click on the Folder icon and find your class or research professor’s folder
- Double-click on the folder with your name.
- Double click on the experiment file
- To get the Folder icon to go back up a step, click on the Folder icon again, then click ONCE only on the little icon that shows an arrow up

16. Getting the last sample out and replacing with a Lock Sample (if auto-eject isn’t turned on)

- a. In “Spectrometer” mode, display “zones” map ()
- b. Right click on sample 48 => select “Sample in Magnet” (3rd choice from the bottom) => OK.

17. Logout: Click “Logout” button underneath spectrum-display