

Chemical Equilibrium in Solution



The activity coefficient in general is approximated by

$$\text{DH Limiting Law; } \log \gamma_i = \frac{-509 z_i^2 \sqrt{I}}{1 + \sqrt{I}}$$

$$K_a \approx \frac{[I_3^-]}{[I_2][I^-]} = K_c$$

Meaning of [X]:


$$[X] = \frac{c_x}{c_0} \text{ where } c_x, c_0 \text{ are concentration of X in mol/L}$$

and concentration of the standard state: 1 mol/L


In effect [X] is the concentration value of X, in mol/l and no unit is involved.

The strategy is to establish the reaction equilibrium in the aqueous solution which is in contact with an organic layer.

This establishes a second (partition) equilibrium system in addition to the reaction equilibrium.



For the chemical equilibrium,



The thermodynamic Equilibrium constant is defined as;


$$K_a = \frac{a_{I_3^-}}{a_{I_2} a_{I^-}} \quad a_i = \gamma_i c_i$$

$$\text{and } K_a = \frac{a_{I_3^-}}{a_{I_2} a_{I^-}} = \frac{\gamma_{I_3^-}}{\gamma_{I_2} \gamma_{I^-}} \frac{[I_3^-]}{[I_2][I^-]} \approx \frac{[I_3^-]}{[I_2][I^-]} = K_c$$


$$K_a \approx K_c$$

implies $\gamma_{I_3^-} = \gamma_{I^-}$; $\gamma_{I_2} = 1$ from the above approximation.


$$K_a \approx \frac{[I_3^-]}{[I_2][I^-]} = K_c$$



To determine the $[I_2]$ one can use the quantitative reaction with the thiosulfate ion.



To calculate the concentrations of the three species in aqueous solution direct reaction method is not feasible as it disturbs the reaction equilibrium. Only the total concentration $[I_2(aq) + I_3^-(aq)]$ is easily determinable, not individual values.



$$k = \frac{I_2(aq)}{I_2(org)}$$


At a given temperature partition coefficient k is a constant for a given equilibrium between the two phases.

If k is known or established, one can find the I_2 concentration in organic layer, and use k to find the $I_2(aq)$ concentration.

Mix		CH ₂ Cl ₂ layer titrate 10.00mL		Aqueous layer titrate 25.00 mL		
Run	CH ₂ Cl ₂ (25 mL) Molarity I ₂	Aqueous Layer (100 mL) Molarity KI	Burette Size, mL	Molarity S ₂ O ₃ ²⁻	Burette Size, mL	Molarity S ₂ O ₃ ²⁻
1	0.080	0.0 (DI)	50	0.1	10	0.01
2	0.040	0.0 (DI)	50	0.1	10	0.01
3	0.020	0.0 (DI)	10	0.1	10	0.01
4	0.080	0.15	10	0.1	50	0.1
5	0.040	0.15	10	0.1	10	0.1
6	0.080	0.03	50	0.1	10	0.1

Determining k

Run	CH ₂ Cl ₂ (25 mL) Molarity I ₂	Aqueous Layer (100 mL) Molarity KI	CH ₂ Cl ₂ layer titrate 10.00mL		Aqueous layer titrate 25.00 mL	
			Vol CH ₂ Cl ₂ mL	[S ₂ O ₃ ²⁻] M	Aq. volume mL	[S ₂ O ₃ ²⁻] M
1	0.080	0.0	10.00	0.10	25.00	0.01
2	0.040	0.0	10.00	0.10	25.00	0.01
3	0.020	0.0	10.00	0.10	25.00	0.01





Run	Vol (mL) 0.1M S ₂ O ₃ ²⁻	[I ₂ org]	CH ₂ Cl ₂ layer titrate 10.00mL		Aqueous layer titrate 25.00 mL	
			Vol (mL) 0.01M S ₂ O ₃ ²⁻	[I ₂ aq]	k	
1						
2						
3						
mean k =						

Determining K_c


Run	CH ₂ Cl ₂ (25 mL) Molarity I ₂	Aqueous Layer (100 mL) Molarity KI	CH ₂ Cl ₂ layer titrate 10.00mL		Aqueous layer titrate 25.00 mL	
			Vol CH ₂ Cl ₂ mL	[S ₂ O ₃ ²⁻] M	Aq. volume mL	[S ₂ O ₃ ²⁻] M
4	0.080	0.15	10.00	0.10	25.00	0.10
5	0.040	0.15	10.00	0.10	25.00	0.10
6	0.080	0.03	10.00	0.10	25.00	0.10

Run	CH ₂ Cl ₂ layer (use -10mL pipette)		Aqueous layer (use 25-mL pipette)				K _c
	Vol (mL) 0.1M S ₂ O ₃ ²⁻	[I ₂ org]	Vol (mL) 0.01M S ₂ O ₃ ²⁻	[I ₂ aq] _{total}	[I ₃ ⁻ aq]	[I ⁻ aq]	
4							
5							
6							


$[I^-(aq)] = [I^-(aq)]_0 - [I_3^-(aq)]$



Titration Reaction



or



@ 'endpoint' $\frac{mol I_2}{mol S_2O_3^{2-}} = \frac{1}{2}$

$$[I_2]V_{I_2} = \frac{1}{2} [S_2O_3^{2-}]V_{S_2O_3^{2-}}$$

$$[I_2] = \frac{1}{2} [S_2O_3^{2-}] \frac{V_{S_2O_3^{2-}}}{V_{I_2}}$$

Titration of the aqueous layer of the reaction mixture:

$\Rightarrow [I_2(aq)] + [I_3^-(aq)]$ of the equilibrium system

Determining K_c

Run	CH ₂ Cl ₂ (25 mL) Molarity I ₂	Aqueous Layer (100 mL) Molarity KI	CH ₂ Cl ₂ layer titrate 10.00mL		Aqueous layer titrate 25.00 mL	
			Vol CH ₂ Cl ₂ mL	[S ₂ O ₃ ²⁻] M	Aq. volume mL	[S ₂ O ₃ ²⁻] M
4	0.080	0.15	10.00	0.10	25.00	0.10
5	0.040	0.15	10.00	0.10	25.00	0.10
6	0.080	0.03	10.00	0.10	25.00	0.10

Run	CH ₂ Cl ₂ layer (use -10mL pipette)		Aqueous layer (use 25-mL pipette)				K _c = $\frac{[I_3^-]_{aq}}{[I_2]_{aq}[I^-]_{aq}}$
	Vol (mL) 0.1M S ₂ O ₃ ²⁻	[I ₂ org]	Vol (mL) 0.01M S ₂ O ₃ ²⁻	[I ₂ aq] _{total}	[I ₂ aq]	[I ₃ ⁻ aq]	
4							
5							
6							

$[I^-(aq)] = [I^-(aq)]_0 - [I_3^-(aq)]$